

Product Name: Recombinant Human SLAMF8 (C-6His)
Catalog #: PHH1539




概述 (Summary)

英文全称	SLAMF8/SLAM family member 8/BLAME
纯度 (Purity)	Greater than 95% as determined by reducing SDS-PAGE
内毒素 (Endotoxin level)	<1 EU/μg as determined by LAL test.
蛋白构建 (Construction)	Recombinant Human SLAM Family Member 8 is produced by our Mammalian expression system and the target gene encoding Ala23-Asp233 is expressed with a 6His tag at the C-terminus.
Accession #	Q9P0V8
蛋白标签 (Tag)	
表达宿主 (Host)	Human Cells
种属 (Species)	Human
预测分子量 (Predicted MW)	26.2 KDa
蛋白形态 (Form)	Lyophilized from a 0.2 μm filtered solution of PBS, pH 7.4.
储存缓冲液 (Buffer)	
运输方式 (Shipping)	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature listed below.
稳定性&储存 (Stability &Storage)	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
复溶 (Reconstitution)	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100μg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles. Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100μg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

电泳图 (SDS-PAGE image)

Product Name: Recombinant Human SLAMF8 (C-6His)
Catalog #: PHH1539



背景 (Background)

分子别名 (Alternative Names)

SLAM family member 8; B-lymphocyte activator macrophage expressed; BCM-like membrane protein; CD353; SLAMF8; BLAME

背景介绍 (References)

SLAM family member 8 (SLAMF8) is a single-pass type I membrane protein and contains 1 Ig-like C2-type domain. SLAMF8 is a member of the CD2 family of cell surface proteins involved in lymphocyte activation. These proteins are characterized by Ig domains and studies of a similar protein in mouse suggest that it may function during B cell lineage commitment. SLAMF8 is expressed in lymph node, spleen, thymus and bone marrow. It may play a role in B-lineage commitment and/or modulation of signaling through the B-cell receptor.

注意事项 (Note)

For Research Use Only , Not for Diagnostic Use.