


Product Name: Recombinant Human ST6GalNAc2 (C-6His)
Catalog #: PHH1563

概述 (Summary)

英文全称	ST6 Sialyltransferase 2/ST6GALNAC2
纯度 (Purity)	Greater than 95% as determined by reducing SDS-PAGE
内毒素 (Endotoxin level)	<1 EU/μg as determined by LAL test.
蛋白构建 (Construction)	Recombinant Human ST6 Sialyltransferase 2 is produced by our Mammalian expression system and the target gene encoding Ser29-Arg374 is expressed with a 6His tag at the C-terminus.
Accession #	Q9UJ37
蛋白标签 (Tag)	
表达宿主 (Host)	Human Cells
种属 (Species)	Human
预测分子量 (Predicted MW)	38.84 KDa
蛋白形态 (Form)	Supplied as a 0.2 μm filtered solution of 20mM PB, 150mM NaCl, 5mM EDTA, 5% Trehalosen, pH 7.4.
储存缓冲液 (Buffer)	
运输方式 (Shipping)	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
稳定性&储存 (Stability &Storage)	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
复溶 (Reconstitution)	

电泳图 (SDS-PAGE image)



背景 (Background)

Product Name: Recombinant Human ST6GalNAc2 (C-6His)
Catalog #: PHH1563



分子别名 (Alternative Names)

Alpha-N-Acetylgalactosaminide Alpha-2; 6-Sialyltransferase 2; GalNAc Alpha-2; 6-Sialyltransferase II; ST6GalNAc II; ST6GalNAcII; SThM; Sialyltransferase 7B; SIAT7-B; ST6GALNAC2; SIAT7B; SIATL1; STHM

背景介绍 (References)

Alpha-N-Acetylgalactosaminide α -2,6-Sialyltransferase 2 (ST6GALNAC2) belongs to the glycosyltransferase 29 family that adds sialic acids to the non-reducing ends of glycoconjugates. At the cell surface, these modifications play roles in cell-cell and cell-substrate interactions, bacterial adhesion and protein targeting. ST6GALNAC2 is localized to the Golgi apparatus membrane. ST6GALNAC2 is highly expressed in lactating mammary glands and the adult testis; it is expressed at lower levels in the kidney. ST6GALNAC2 can catalyze 2,6-sialylation of the Tn antigen.

注意事项 (Note)

For Research Use Only , Not for Diagnostic Use.