


Product Name: Recombinant Human SUMF1 (C-6His)
Catalog #: PHH1579



概述 (Summary)

英文全称	Sulfatase Modifying Factor 1/SUMF1
纯度 (Purity)	Greater than 95% as determined by reducing SDS-PAGE
内毒素 (Endotoxin level)	<1 EU/μg as determined by LAL test.
蛋白构建 (Construction)	Recombinant Human Sulfatase Modifying Factor 1 is produced by our Mammalian expression system and the target gene encoding Ser34-Asp374 is expressed with a 6His tag at the C-terminus.
Accession #	Q8NBK3
蛋白标签 (Tag)	
表达宿主 (Host)	Human Cells
种属 (Species)	Human
预测分子量 (Predicted MW)	38.4 KDa
蛋白形态 (Form)	Supplied as a 0.2 μm filtered solution of 20mM Tris-HCl, 150mM NaCl, 2mM CaCl ₂ , 10% Glycerol, pH 7.5.
储存缓冲液 (Buffer)	
运输方式 (Shipping)	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
稳定性&储存 (Stability &Storage)	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
复溶 (Reconstitution)	

电泳图 (SDS-PAGE image)



背景 (Background)

Product Name: Recombinant Human SUMF1 (C-6His)
Catalog #: PHH1579



分子别名 (Alternative Names)

Sulfatase-Modifying Factor 1; C-Alpha-Formylglycine-Generating Enzyme 1; SUMF1; FGE

背景介绍 (References)

Human Sulfatase Modifying Factor 1 (SUMF1) is a 42kDa protein. SUMF1 is a Ca²⁺-binding member of the sulfatase-modifying factor family. SUMF1 is a soluble ER luminal glycoprotein, it converts inactive sulfatases into an active form by transforming a catalytic site cysteine into a formylglycine residue. In the ER, SUMF1 can exist as either a monomer, or a disulfide-linked homodimer or a heterodimer with SUMF2. Three splice isoforms are known.

注意事项 (Note)

For Research Use Only , Not for Diagnostic Use.