


## 概述 (Summary)

|                             |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
|-----------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 英文全称                        | TFF2/Trefoil factor 2                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| 纯度 (Purity)                 | Greater than 95% as determined by reducing SDS-PAGE                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| 内毒素 (Endotoxin level)       | <1 EU/μg as determined by LAL test.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| 蛋白构建 (Construction)         | Recombinant Human Trefoil Factor 2 is produced by our Mammalian expression system and the target gene encoding Glu24-Tyr129 is expressed with a 6His tag at the C-terminus.                                                                                                                                                                                                                                                                                                                                                                                                   |
| Accession #                 | Q03403                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| 蛋白标签 (Tag)                  |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| 表达宿主 (Host)                 | Human Cells                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| 种属 (Species)                | Human                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| 预测分子量 (Predicted MW)        | 13 KDa                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| 蛋白形态 (Form)                 | Lyophilized from a 0.2 μm filtered solution of PBS, pH 7.4.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| 储存缓冲液 (Buffer)              |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| 运输方式 (Shipping)             | The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature listed below.                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| 稳定性&储存 (Stability &Storage) | Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.                                                                                                                                                                                                                                                                                                                                                                                                          |
| 复溶 (Reconstitution)         | Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100μg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles. Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100μg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles. |

## 电泳图 (SDS-PAGE image)

**Product Name: Recombinant Human TFF2 (C-6His)**  
**Catalog #: PHH1618**



## 背景 (Background)

### 分子别名 (Alternative Names)

Trefoil Factor 2; Spasmolysin; Spasmolytic Polypeptide; SP; TFF2; SML1

### 背景介绍 (References)

Trefoil Factor 2 (TFF2) is a member of the trefoil family and contains two P-type (trefoil) domains. Members of this family are characterized by having at least one copy of the trefoil motif, a 40-amino acid domain that contains three conserved disulfides. TFF2 is a secreted protein and specifically expressed in the stomach. TFF2 inhibits gastrointestinal motility and gastric acid secretion. TFF2 could function as a structural component of gastric mucus, possibly by stabilizing glycoproteins in the mucus gel through interactions with carbohydrate side chains.

## 注意事项 (Note)

For Research Use Only , Not for Diagnostic Use.