


产品名称: Recombinant Human TINAGL1 (C-6His)  
产品货号: PHH1657



## 概述 (Summary)

英文全称	TINAGL1/Tubulointerstitial nephritis antigen-like
纯度 (Purity)	Greater than 95% as determined by reducing SDS-PAGE
内毒素 (Endotoxin level)	<1 EU/ $\mu$ g as determined by LAL test.
蛋白构建 (Construction)	Recombinant Human Tubulointerstitial Nephritis Antigen-like 1 is produced by our Mammalian expression system and the target gene encoding Ala22-His467 is expressed with a 6His tag at the C-terminus.
Accession #	Q9GZM7
表达宿主 (Host)	Human Cells
种属 (Species)	Human
预测分子量 (Predicted MW)	51.1 KDa
制剂 (Form)	Supplied as a 0.2 $\mu$ m filtered solution of 20mM Tris-HCl, 150mM NaCl, 10% Glycerol, 0.05% Brij35, pH 7.5.
运输方式 (Shipping)	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
稳定性&储存 (Stability &Storage)	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
复溶 (Reconstitution)	

## 电泳图 (SDS-PAGE image)



## 背景 (Background)

产品名称: Recombinant Human TINAGL1 (C-6His)  
产品货号: PHH1657



---

**分子别名 (Alternative Names)**

Tubulointerstitial nephritis antigen-like (TINAGL1); Glucocorticoid-inducible protein 5; Oxidized LDL-responsive gene 2 protein; TIN Ag-related protein; Tubulointerstitial nephritis antigen-related protein; GIS5; LCN7; OLRG2 and TINAGL

**背景介绍 (References)**

TINAGL1 is a secreted protein and contains one SMB (somatomedin-B) domain. TINAGL1 is expressed in aorta, lymph nodes, heart, bone marrow, placenta, lung, kidney, skeletal muscle, pancreas, thyroid, adrenal gland and so on. TINAGL1 may be implicated in the adrenocortical zonation and in mechanisms for repressing the CYP11B1 gene expression in adrenocortical cells.

**注意事项 (Note)**

For research use only .