


Product Name: Recombinant Human HATL5 (C-6His)
Catalog #: PHH1672



概述 (Summary)

英文全称	TMPRSS11B/Transmembrane protease serine 11B
纯度 (Purity)	Greater than 95% as determined by reducing SDS-PAGE
内毒素 (Endotoxin level)	<1 EU/μg as determined by LAL test.
蛋白构建 (Construction)	Recombinant Human Airway Trypsin-Like Protease 5 is produced by our Mammalian expression system and the target gene encoding Leu39-Leu416 is expressed with a 6His tag at the C-terminus.
Accession #	Q86T26
蛋白标签 (Tag)	
表达宿主 (Host)	Human Cells
种属 (Species)	Human
预测分子量 (Predicted MW)	43.1 KDa
蛋白形态 (Form)	Supplied as a 0.2 μm filtered solution of 20mM PB, 150mM NaCl, pH 7.4.
储存缓冲液 (Buffer)	
运输方式 (Shipping)	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
稳定性&储存 (Stability &Storage)	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
复溶 (Reconstitution)	

电泳图 (SDS-PAGE image)



背景 (Background)

分子别名 (Alternative Names) Transmembrane Protease Serine 11B; Airway Trypsin-Like Protease 5;

Product Name: Recombinant Human HATL5 (C-6His)
Catalog #: PHH1672



背景介绍 (References)

TMPRSS11B; HATL5

Transmembrane Protease Serine 11B (TMPRSS11B) is a single-pass type II membrane protein member of the peptidase S1 family. TMPRSS11B contains one peptidase S1 domain and one SEA domain. TMPRSS11B is a serine protease that may play some biological role in the host defense system on the mucous membrane independently of or in cooperation with other substances in airway mucous or bronchial secretions.

注意事项 (Note)

For Research Use Only , Not for Diagnostic Use.