


概述 (Summary)

英文全称	VSIG8/V-Set and Ig Domain-Containing Protein 8
纯度 (Purity)	Greater than 95% as determined by reducing SDS-PAGE
内毒素 (Endotoxin level)	<1 EU/μg as determined by LAL test.
蛋白构建 (Construction)	Recombinant Human V-Set And Ig Domain-Containing Protein 8 is produced by our Mammalian expression system and the target gene encoding Val22-Gly263 is expressed with a 6His tag at the C-terminus.
Accession #	P0DPA2
蛋白标签 (Tag)	
表达宿主 (Host)	Human Cells
种属 (Species)	Human
预测分子量 (Predicted MW)	27.9 KDa
蛋白形态 (Form)	Lyophilized from a 0.2 μm filtered solution of 20mM NaAC-HAc, 5%Trehalose, pH 4.5.
储存缓冲液 (Buffer)	
运输方式 (Shipping)	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature listed below.
稳定性&储存 (Stability &Storage)	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
复溶 (Reconstitution)	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100μg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles. Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100μg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

电泳图 (SDS-PAGE image)

Product Name: Recombinant Human VSIG8 (C-6His)
Catalog #: PHH1820



背景 (Background)

分子别名 (Alternative Names)

V-set and immunoglobulin domain-containing protein 8; VSIG8; C1orf204

背景介绍 (References)

V-set and immunoglobulin domain-containing protein 8(VSIG8) is a single-pass type I membrane protein. The human VSIG8 cDNA encodes 414 amino acids (aa) including a 21 aa signal sequence, a 242 aa extracellular domain (ECD) containing 2 Ig-like V-type (immunoglobulin-like) domains, a 21 aa transmembrane domain and a 130 aa cytoplasmic domain. The function of VSIG8 is not clear.

注意事项 (Note)

For Research Use Only , Not for Diagnostic Use.