


概述 (Summary)

英文全称	VSIG8/V-Set and Ig Domain-Containing Protein 8
纯度 (Purity)	Greater than 95% as determined by reducing SDS-PAGE
内毒素 (Endotoxin level)	<1 EU/μg as determined by LAL test.
蛋白构建 (Construction)	Recombinant Human V-Set And Ig Domain-Containing Protein 8 is produced by our Mammalian expression system and the target gene encoding Val22-Gly263 is expressed with a human IgG1 Fc tag at the C-terminus.
Accession #	P0DPA2
蛋白标签 (Tag)	
表达宿主 (Host)	Human Cells
种属 (Species)	Human
预测分子量 (Predicted MW)	54.2 KDa
蛋白形态 (Form)	Lyophilized from a 0.2 μm filtered solution of 20mM PB, 150mM NaCl, pH 7.4.
储存缓冲液 (Buffer)	
运输方式 (Shipping)	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature listed below.
稳定性&储存 (Stability &Storage)	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
复溶 (Reconstitution)	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100μg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles. Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100μg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

电泳图 (SDS-PAGE image)

Product Name: Recombinant Human VSIG8 (C-Fc)
Catalog #: PHH1822



背景 (Background)

分子别名 (Alternative Names)

V-set and immunoglobulin domain-containing protein 8;VSIG8;C1orf204

背景介绍 (References)

V-set and immunoglobulin domain-containing protein 8(VSIG8) is a single-pass type I membrane protein. The human VSIG8 cDNA encodes 414 amino acids (aa) including a 21 aa signal sequence, a 242 aa extracellular domain (ECD) containing 2 Ig-like V-type (immunoglobulin-like) domains, a 21 aa transmembrane domain and a 130 aa cytoplasmic domain. The function of VSIG8 is not clear.

注意事项 (Note)

For Research Use Only , Not for Diagnostic Use.