


概述 (Summary)

英文全称	VSTM2A/V-set and transmembrane domain-containing protein 2A
纯度 (Purity)	Greater than 95% as determined by reducing SDS-PAGE
内毒素 (Endotoxin level)	<1 EU/μg as determined by LAL test.
蛋白构建 (Construction)	Recombinant Human V-set And Transmembrane Domain-Containing Protein 2A is produced by our Mammalian expression system and the target gene encoding Ser24-Phe243 is expressed with a 6His tag at the C-terminus.
Accession #	AAH28404.1
蛋白标签 (Tag)	
表达宿主 (Host)	Human Cells
种属 (Species)	Human
预测分子量 (Predicted MW)	25.5 KDa
蛋白形态 (Form)	Lyophilized from a 0.2 μm filtered solution of 20mM PB, 150mM NaCl, pH 7.4.
储存缓冲液 (Buffer)	
运输方式 (Shipping)	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature listed below.
稳定性&储存 (Stability &Storage)	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
复溶 (Reconstitution)	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100μg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles. Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100μg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

电泳图 (SDS-PAGE image)

Product Name: Recombinant Human VSTM2A (C-6His)
Catalog #: PHH1828



背景 (Background)

分子别名 (Alternative Names)

V-set and transmembrane domain-containing protein 2A;VSTM2

背景介绍 (References)

VSTM2A, also called VSTM2, is short for V-set and transmembrane domain-containing protein 2A. It is a 236 aa. protein which can be departed into two parts. The first 24 aa. is a signal peptide, and the last is a V-set and transmembrane domain-containing protein 2A. There is a natural variant location which is E84K.

注意事项 (Note)

For Research Use Only , Not for Diagnostic Use.