

**Product Name: Recombinant Human IL-2 Superkine (Mutant, C-6His)**  
**Catalog #: PHH1918**

---


**概述 (Summary)**

英文全称	IL-2/Interleukin-2/T cell growth factor/TCGF (Mutant)
纯度 (Purity)	Greater than 95% as determined by reducing SDS-PAGE
内毒素 (Endotoxin level)	<1 EU/μg as determined by LAL test.
蛋白构建 (Construction)	Recombinant Human Interleukin-2 Superkine is produced by our Mammalian expression system and the target gene encoding Ala21-Thr153 (L100F, R101D, L105V, I106V, I112F) is expressed with a 6His tag at the C-terminus.
Accession #	P60568
蛋白标签 (Tag)	
表达宿主 (Host)	Human Cells
种属 (Species)	Human
预测分子量 (Predicted MW)	16.5 KDa
蛋白形态 (Form)	Lyophilized from a 0.2 μm filtered solution of 4mM HCl.
储存缓冲液 (Buffer)	
运输方式 (Shipping)	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature listed below.
稳定性&储存 (Stability &Storage)	Lyophilized protein should be stored at ≤ -20°C, stable for one year after receipt. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at ≤ -20°C for 3 months.
复溶 (Reconstitution)	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100μg/ml. Dissolve the lyophilized protein in 4mM HCl. Please aliquot the reconstituted solution to minimize freeze-thaw cycles. Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100μg/ml. Dissolve the lyophilized protein in 4mM HCl. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

**电泳图 (SDS-PAGE image)**

**Product Name: Recombinant Human IL-2 Superkine (Mutant, C-6 His)**  
**Catalog #: PHH1918**

---



## 背景 (Background)

### 分子别名 (Alternative Names)

Interleukin-2; IL-2; T-Cell Growth Factor; TCGF; Aldesleukin; IL2

### 背景介绍 (References)

Interleukin-2(IL-2) is an interleukin, a type of cytokine signaling molecule in the immune system , belongs to the IL-2 family. It is a powerful immunoregulatory lymphokine produced by T-cells in response to antigenic or mitogenic stimulation. IL-2/IL-2R signaling is required for T-cell proliferation and other fundamental functions that are essential for the immune response. IL-2 stimulates growth and differentiation of B-cells, NK cells, lymphokine-activated killer cells, monocytes, macrophages and oligodendrocytes. New research has shown that IL-2 mutant reduced toxicity while being more potent at stimulating anti-tumor effector immune cells.

## 注意事项 (Note)

For Research Use Only , Not for Diagnostic Use.