


概述 (Summary)

英文全称	REG3G/Regenerating islet-derived protein 3-gamma/ Reg3G/Reg III-gamma
纯度 (Purity)	Greater than 95% as determined by reducing SDS-PAGE
内毒素 (Endotoxin level)	<1 EU/μg as determined by LAL test.
蛋白构建 (Construction)	Recombinant Human Regenerating Islet-derived Protein 3-gamma is produced by our Mammalian expression system and the target gene encoding Glu27-Asp175 is expressed
Accession #	Q6UW15
蛋白标签 (Tag)	
表达宿主 (Host)	Human Cells
种属 (Species)	Human
预测分子量 (Predicted MW)	46.5 KDa
蛋白形态 (Form)	Supplied as a 0.2 μm filtered solution of 20 mM PB, 500 mM NaCl, 10% Glycerol, 0.05% Tween 80, pH7.5.
储存缓冲液 (Buffer)	
运输方式 (Shipping)	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
稳定性&储存 (Stability &Storage)	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
复溶 (Reconstitution)	

电泳图 (SDS-PAGE image)



背景 (Background)

Product Name: Recombinant Human REG3G (C-Fc)
Catalog #: PHH1985



分子别名 (Alternative Names)

Regenerating islet-derived protein 3-gamma; REG-3-gamma; Pancreatitis-associated protein 1B; REG3G; PAP-1B; Regenerating islet-derived protein III-gamma

背景介绍 (References)

Pancreatitis-associated protein 1B, Regenerating islet-derived protein III-gamma, is a secreted protein which contains 1 C-type lectin domain. It is expressed almost exclusively in the pancreas and also expressed in testis, but not found in small intestine. This protein might be a stress protein involved in the control of bacterial proliferation.

注意事项 (Note)

For Research Use Only , Not for Diagnostic Use.