


## 概述 (Summary)

|                             |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
|-----------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 英文全称                        | COL21A1/Collagen alpha-1(XXI) chain                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| 纯度 (Purity)                 | Greater than 95% as determined by reducing SDS-PAGE                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| 内毒素 (Endotoxin level)       | <1 EU/μg as determined by LAL test.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| 蛋白构建 (Construction)         | Recombinant Human Collagen Alpha-1(XXI) Chain is produced by our Mammalian expression system and the target gene encoding Glu23-Tyr957 is expressed with a 6His tag at the C-terminus.                                                                                                                                                                                                                                                                                                                                                                                        |
| Accession #                 | Q96P44                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| 蛋白标签 (Tag)                  |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| 表达宿主 (Host)                 | Human Cells                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| 种属 (Species)                | Human                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| 预测分子量 (Predicted MW)        | 97.9 KDa                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| 蛋白形态 (Form)                 | Lyophilized from a 0.2 μm filtered solution of PBS, pH 7.4.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| 储存缓冲液 (Buffer)              |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| 运输方式 (Shipping)             | The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature listed below.                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| 稳定性&储存 (Stability &Storage) | Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.                                                                                                                                                                                                                                                                                                                                                                                                          |
| 复溶 (Reconstitution)         | Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100μg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles. Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100μg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles. |

## 电泳图 (SDS-PAGE image)

**Product Name: Recombinant Human COL21A1 (C-6His)**  
**Catalog #: PHH2007**



## 背景 (Background)

### 分子别名 (Alternative Names)

COL21A1; collagen, type XXI, alpha 1; collagen alpha-1(XXI) chain; alpha 1 type XXI collagen; FP633; COLA1L; dJ708F5.1; dJ682J15.1; FLJ39125; FLJ44623; MGC26619; DKFZp564B052

### 背景介绍 (References)

Collagen  $\alpha$ -1(XXI) Chain (COL21A1) is a member of the fibril-associated collagens with interrupted helices (FACIT) family. COL21A1 is a secreted protein and contains six collagen-like domains, one TSP N-terminal (TSPN) domain, and one VWFA domain. COL21A1 is widely expressed in many tissues with the highest expression observed at the fetal stage. COL21A1 is stimulated by PDGF/platelet-derived growth factor. Type XXI collagen is localized to tissues containing type I collagen; it may serve to maintain the integrity of the extracellular matrix.

## 注意事项 (Note)

For Research Use Only , Not for Diagnostic Use.