


Product Name: Recombinant Human Carbonic Anhydrase 9 (C-6His)
Catalog #: PHH2186

概述 (Summary)

英文全称	Carbonic Anhydrase IX/CA9
纯度 (Purity)	Greater than 95% as determined by reducing SDS-PAGE
内毒素 (Endotoxin level)	<1 EU/μg as determined by LAL test.
蛋白构建 (Construction)	Recombinant Human Carbonic Anhydrase 9 is produced by our Mammalian expression system and the target gene encoding Gln38-Asp414 is expressed with a 6His tag at the C-terminus.
Accession #	Q16790
蛋白标签 (Tag)	
表达宿主 (Host)	Human Cells
种属 (Species)	Human
预测分子量 (Predicted MW)	42 KDa
蛋白形态 (Form)	Lyophilized from a 0.2 μm filtered solution of PBS, pH 7.4.
储存缓冲液 (Buffer)	
运输方式 (Shipping)	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature listed below.
稳定性&储存 (Stability &Storage)	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
复溶 (Reconstitution)	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100μg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles. Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100μg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

电泳图 (SDS-PAGE image)

Product Name: Recombinant Human Carbonic Anhydrase 9 (C-6His)
Catalog #: PHH2186



背景 (Background)

分子别名 (Alternative Names)

CA9; CA-IX; Carbonic Anhydrase IX; Carbonate dehydratase IX; G250; MN; P54/58N; RCC; RCC-associated protein G250

背景介绍 (References)

Carbonic anhydrases IX (CA IX), also known as membrane antigen MN or CA9, is a member of the carbonic anhydrase (CA) family and may be involved in cell proliferation and cellular transformation. CAs are zinc metalloenzymes that catalyze the reversible hydration of carbon dioxide ($\text{H}_2\text{O} + \text{CO}_2 = \text{H}^+ + \text{HCO}_3^-$) and thus participate in a variety of biological and physical processes. CA9 is a transmembrane enzyme expressed primarily in carcinoma cells. It is one of the best markers for hypoxia and for RCC. Appears to be a novel specific biomarker for a cervical neoplasia.

注意事项 (Note)

For Research Use Only , Not for Diagnostic Use.