


## 概述 (Summary)

|                             |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
|-----------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 英文全称                        | HER3/Receptor Tyrosine-Protein Kinase ErbB-3/ERBB3 (Ser20-Thr643)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| 纯度 (Purity)                 | Greater than 95% as determined by reducing SDS-PAGE                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| 内毒素 (Endotoxin level)       | <1 EU/μg as determined by LAL test.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| 蛋白构建 (Construction)         | Recombinant Human Receptor Tyrosine-Protein Kinase ErbB-3 is produced by our Mammalian expression system and the target gene encoding Ser20-Thr643 is expressed with a human IgG1 Fc tag at the C-terminus.                                                                                                                                                                                                                                                                                                                                                                   |
| Accession #                 | P21860                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| 蛋白标签 (Tag)                  |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| 表达宿主 (Host)                 | Human Cells                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| 种属 (Species)                | Human                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| 预测分子量 (Predicted MW)        | 95.6 KDa                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| 蛋白形态 (Form)                 | Lyophilized from a 0.2 μm filtered solution of 20 mM Tris-HCl, 10% Trehalose, 2% Mannitol, 0.05% Tween 80, pH 9.0.                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| 储存缓冲液 (Buffer)              |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| 运输方式 (Shipping)             | The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature listed below.                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| 稳定性&储存 (Stability &Storage) | Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.                                                                                                                                                                                                                                                                                                                                                                                                          |
| 复溶 (Reconstitution)         | Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100μg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles. Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100μg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles. |

## 电泳图 (SDS-PAGE image)

**Product Name: Recombinant Human HER3 (C-Fc)**  
**Catalog #: PHH2269**



## 背景 (Background)

### 分子别名 (Alternative Names)

Proto-oncogene-like protein c-ErbB-3; Tyrosine kinase-type cell surface receptor HER3; ERBB3; HER3

### 背景介绍 (References)

Receptor tyrosine-protein kinase erbB-3 is an enzyme that in humans is encoded by the ERBB3 gene. This gene encodes a member of the epidermal growth factor receptor (EGFR) family of receptor tyrosine kinases. ERBB3 belongs to the protein kinase superfamily, tyr protein kinase family and EGF receptor subfamily. It contains 1 protein kinase domain and it is expressed in Epithelial tissues and brain. This membrane-bound protein has a neuregulin binding domain but not an active kinase domain. It therefore can bind this ligand but not convey the signal into the cell through protein phosphorylation. However, it does form heterodimers with other EGF receptor family members which do have kinase activity.

## 注意事项 (Note)

For Research Use Only , Not for Diagnostic Use.