


Product Name: Recombinant Human CD39L1 (C-6His)
Catalog #: PHH2374



概述 (Summary)

英文全称	CD39L1/Ecto-Nucleoside Triphosphate Diphosphohydrolase 2
纯度 (Purity)	Greater than 95% as determined by reducing SDS-PAGE
内毒素 (Endotoxin level)	<1 EU/μg as determined by LAL test.
蛋白构建 (Construction)	Recombinant Human Ecto-Nucleoside Triphosphate Diphosphohydrolase 2 is produced by our Mammalian expression system and the target gene encoding Thr29-Asp460 is expressed with a 6His tag at the C-terminus.
Accession #	Q9Y5L3
蛋白标签 (Tag)	
表达宿主 (Host)	Human Cells
种属 (Species)	Human
预测分子量 (Predicted MW)	48 KDa
蛋白形态 (Form)	Supplied as a 0.2 μm filtered solution of 20mM Tris-HCl, 150mM NaCl, 5mM CaCl ₂ , 10% Glycerol, pH 7.4.
储存缓冲液 (Buffer)	
运输方式 (Shipping)	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
稳定性&储存 (Stability &Storage)	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
复溶 (Reconstitution)	

电泳图 (SDS-PAGE image)



背景 (Background)

Product Name: Recombinant Human CD39L1 (C-6His)
Catalog #: PHH2374



分子别名 (Alternative Names)

Ectonucleoside triphosphate diphosphohydrolase 2; Entpd2

背景介绍 (References)

CD39L1 protein (ENTPD2 or NTPDase2) is a member of the ecto-nucleoside triphosphate diphosphohydrolase family which the main role is termination of purinergic signaling. CD39L1 gene encodes a precursor protein with 495 amino acid residues which generates a 437 amino acid residues mature protein after processing. It is an ecto-nucleotidase that found on the surface of vascular adventitial cells and accessory vascular cells. CD39L1 is a Ca^{2+} - and Mg^{2+} -dependent enzyme that activates platelets by preferentially converting ATP to ADP. CD39L1 plays a role in regulating thrombosis and inflammation which is considered to be a therapeutic target for thromboregulation and the treatment of vascular inflammation. Alternative splicing of CD39L1 gene results in multiple transcript variants.

注意事项 (Note)

For Research Use Only , Not for Diagnostic Use.