


概述 (Summary)

英文全称	PSMA/FOLH1/GCP2/FGCP/GCP11/mGCP/NAALADase I
纯度 (Purity)	Greater than 95% as determined by reducing SDS-PAGE
内毒素 (Endotoxin level)	Please contact with the lab for this information
蛋白构建 (Construction)	Recombinant Human Glutamate carboxypeptidase 2 is produced by our Mammalian expression system and the target gene encoding Lys44-Ala750 is expressed with a human IgG1 Fc tag at the N-terminus.
Accession #	AAC83972.1
蛋白标签 (Tag)	
表达宿主 (Host)	Human Cells
种属 (Species)	Human
预测分子量 (Predicted MW)	105.3 kDa
蛋白形态 (Form)	Lyophilized protein should be stored at $\leq -20^{\circ}\text{C}$, stable for one year after receipt. Reconstituted protein solution can be stored at $2-8^{\circ}\text{C}$ for 2-7 days. Aliquots of reconstituted samples are stable at $\leq -20^{\circ}\text{C}$ for 3 months.
储存缓冲液 (Buffer)	
运输方式 (Shipping)	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
稳定性&储存 (Stability &Storage)	0.0
复溶 (Reconstitution)	

电泳图 (SDS-PAGE image)



背景 (Background)

Product Name: Recombinant Human PSMA (N-Fc)
Catalog #: PHH2436



分子别名 (Alternative Names)	0.0
背景介绍 (References)	0.0

注意事项 (Note)

For Research Use Only , Not for Diagnostic Use.