


概述 (Summary)

英文全称	EDIL3/EGF-like repeats and discoidin I-like domains 3/Del-1/Integrin-binding protein DEL1
纯度 (Purity)	Greater than 95% as determined by reducing SDS-PAGE
内毒素 (Endotoxin level)	Please contact with the lab for this information
蛋白构建 (Construction)	Recombinant Human EGF-Like Repeat and Discoidin I-Like Domain-Containing Protein 3 is produced by our Mammalian expression system and the target gene encoding Asp24-Glu480 is expressed with a human IgG1 Fc tag at the C-terminus.
Accession #	O43854-1
蛋白标签 (Tag)	
表达宿主 (Host)	Human Cells
种属 (Species)	Human
预测分子量 (Predicted MW)	78.3 kDa
蛋白形态 (Form)	Lyophilized protein should be stored at $\leq -20^{\circ}\text{C}$, stable for one year after receipt. Reconstituted protein solution can be stored at $2-8^{\circ}\text{C}$ for 2-7 days. Aliquots of reconstituted samples are stable at $\leq -20^{\circ}\text{C}$ for 3 months.
储存缓冲液 (Buffer)	
运输方式 (Shipping)	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
稳定性&储存 (Stability &Storage)	0.0
复溶 (Reconstitution)	

电泳图 (SDS-PAGE image)



Product Name: Recombinant Human EDIL3 (C-Fc)
Catalog #: PHH2439



背景 (Background)

分子别名 (Alternative Names)	0.0
背景介绍 (References)	0.0

注意事项 (Note)

For Research Use Only , Not for Diagnostic Use.