


Product Name: Recombinant Human DLL3 (C-6His)
Catalog #: PHH2462



概述 (Summary)

英文全称	DLL3
纯度 (Purity)	Greater than 95% as determined by reducing SDS-PAGE
内毒素 (Endotoxin level)	<1 EU/μg as determined by LAL test.
蛋白构建 (Construction)	Recombinant Human Delta-like protein 3 is produced by our Mammalian expression system and the target gene encoding Ala27-Leu492 is expressed with a 6His tag at the C-terminus.
Accession #	Q9NYJ7
蛋白标签 (Tag)	
表达宿主 (Host)	Human Cells
种属 (Species)	Human
预测分子量 (Predicted MW)	49.4 KDa
蛋白形态 (Form)	Supplied as a 0.2 μm filtered solution of 50 mM Tris-HCl, 150 mM NaCl, 0.1 M Glycine, 0.2 M Arg, 5% Trehalose, pH7.5.
储存缓冲液 (Buffer)	
运输方式 (Shipping)	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
稳定性&储存 (Stability &Storage)	Store at ≤ -70°C, stable for 6 months after receipt. Store at ≤ -70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
复溶 (Reconstitution)	

电泳图 (SDS-PAGE image)



背景 (Background)

Product Name: Recombinant Human DLL3 (C-6His)
Catalog #: PHH2462



分子别名 (Alternative Names)

Delta-like protein 3; Drosophila Delta homolog 3; Delta3; DLL3

背景介绍 (References)

Delta-like protein 3 (DLL3) is a transmembrane protein that belongs to the Delta/Serrate/Lag-2 (DSL) family of Notch ligands. DLL3 is highly expressed in tumor cells and could be involved in the modification of tumor (immune) microenvironment. Neuroendocrine carcinomas are very aggressive tumors with few treatment options. DLL3 seems to be an optimal target for therapeutic intervention, as it is expressed mainly on the membrane of tumor cells with neuroendocrine origin.

注意事项 (Note)

For Research Use Only , Not for Diagnostic Use.