


## 概述 (Summary)

|                             |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
|-----------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 英文全称                        | Fc gamma RIIB/CD32b/Low affinity immunoglobulin gamma Fc region receptor II                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| 纯度 (Purity)                 | Greater than 95% as determined by reducing SDS-PAGE                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| 内毒素 (Endotoxin level)       | <1 EU/μg as determined by LAL test.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| 蛋白构建 (Construction)         | Recombinant Mouse Low Affinity Immunoglobulin Gamma Fc Region Receptor II is produced by our Mammalian expression system and the target gene encoding Thr30-Pro210 is expressed with a 6His tag at the C-terminus.                                                                                                                                                                                                                                                                                                                                                            |
| Accession #                 | P08101                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| 蛋白标签 (Tag)                  |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| 表达宿主 (Host)                 | Human Cells                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| 种属 (Species)                | Mouse                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| 预测分子量 (Predicted MW)        | 21.6 KDa                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| 蛋白形态 (Form)                 | Lyophilized from a 0.2 μm filtered solution of PBS, pH 7.4.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| 储存缓冲液 (Buffer)              |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| 运输方式 (Shipping)             | The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature listed below.                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| 稳定性&储存 (Stability &Storage) | Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.                                                                                                                                                                                                                                                                                                                                                                                                          |
| 复溶 (Reconstitution)         | Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100μg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles. Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100μg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles. |

## 电泳图 (SDS-PAGE image)

**Product Name: Recombinant Mouse CD32b (C-6His)**  
**Catalog #: PHM1098**



## 背景 (Background)

### 分子别名 (Alternative Names)

Low affinity immunoglobulin gamma Fc region receptor II; Fc gamma receptor IIB; Fc-gamma RII; Fc-gamma-RIIB; FcRII; IgG Fc receptor II beta; Lymphocyte antigen 17; Ly-17; CD32; Fcgr2; Fcgr2b; Ly-17

### 背景介绍 (References)

Low affinity immunoglobulin gamma Fc region receptor II (CD32B) is a single-pass type I membrane protein and contains 2 Ig-like C2-type (immunoglobulin-like) domains. The inhibitory CD32B is expressed on B cells and myeloid dendritic cells. Ligation of CD32B on B cells downregulates antibody production and may, in some circumstances, promote apoptosis. Co-ligation of CD32B on dendritic cells inhibits maturation and blocks cell activation. CD32B may also be a target for monoclonal antibody therapy for malignancies.

## 注意事项 (Note)

For Research Use Only , Not for Diagnostic Use.