


Product Name: Recombinant Mouse 4-1BB (C-Fc)
Catalog #: PHM1696



概述 (Summary)

英文全称	4-1BB/CD137/TNFRSF9/ILA
纯度 (Purity)	Greater than 95% as determined by reducing SDS-PAGE
内毒素 (Endotoxin level)	<1 EU/μg as determined by LAL test.
蛋白构建 (Construction)	Recombinant Mouse 4-1BB ligand receptor is produced by our Mammalian expression system and the target gene encoding Val24-Leu187 is expressed with a human IgG1 Fc tag at the C-terminus.
Accession #	P20334
蛋白标签 (Tag)	
表达宿主 (Host)	Human Cells
种属 (Species)	Mouse
预测分子量 (Predicted MW)	44.7 KDa
蛋白形态 (Form)	Supplied as a 0.2 μm filtered solution of 20mM PB, 150mM NaCl, pH 7.4.
储存缓冲液 (Buffer)	
运输方式 (Shipping)	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
稳定性&储存 (Stability &Storage)	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
复溶 (Reconstitution)	

电泳图 (SDS-PAGE image)



背景 (Background)

分子别名 (Alternative Names)	Secreted CD137 antigen ;Tumor necrosis factor receptor superfamily member
--------------------------	---

Product Name: Recombinant Mouse 4-1BB (C-Fc)
Catalog #: PHM1696



背景介绍 (References)

9 ;Tnfrsf9

Tumor necrosis factor receptor superfamily member 9 (TNFRSF9) is a member of the tumor necrosis factor (TNF) receptor family. It can be induced by lymphocyte activation (ILA) and is expressed by activated T cells, but to a larger extent on CD8 than on CD4 T cells. In addition, TNFRSF9 expression is found on dendritic cells, follicular dendritic cells, natural killer cells, granulocytes and cells of blood vessel walls at sites of inflammation. As receptor for TNFSF9/4-1BBL, it can activate T cells and the cross-linking of this protein can enhance T cell proliferation, IL-2 secretion survival and cytolytic activity. Further, it can enhance immune activity to eliminate tumors in mice.

注意事项 (Note)

For Research Use Only , Not for Diagnostic Use.