


概述 (Summary)

英文全称	VSIG8/V-Set and Ig Domain-Containing Protein 8
纯度 (Purity)	Greater than 95% as determined by reducing SDS-PAGE
内毒素 (Endotoxin level)	<1 EU/μg as determined by LAL test.
蛋白构建 (Construction)	Recombinant Mouse V-set And Immunoglobulin Domain-containing Protein 8 is produced by our Mammalian expression system and the target gene encoding Val22-Gly262 is expressed with a 6His tag at the C-terminus.
Accession #	Q6P3A4
蛋白标签 (Tag)	
表达宿主 (Host)	Human Cells
种属 (Species)	Mouse
预测分子量 (Predicted MW)	27.6 KDa
蛋白形态 (Form)	Supplied as a 0.2 μm filtered solution of PBS, 10% Glycerol, pH 7.4.
储存缓冲液 (Buffer)	
运输方式 (Shipping)	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
稳定性&储存 (Stability &Storage)	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
复溶 (Reconstitution)	

电泳图 (SDS-PAGE image)



背景 (Background)

Product Name: Recombinant Mouse VSIG8 (C-6His)
Catalog #: PHM1821



分子别名 (Alternative Names)

V-set and immunoglobulin domain-containing protein 8; Vsig8

背景介绍 (References)

V-set and immunoglobulin domain-containing protein 8(VSIG8) is a single-pass type I membrane protein. The mouse VSIG8 cDNA encodes 417 amino acids (aa) including a 21 aa signal sequence, a 241 aa extracellular domain (ECD) containing 2 Ig-like V-type (immunoglobulin-like) domains, a 21 aa transmembrane domain and a 134 aa cytoplasmic domain. The function of VSIG8 is not clear.

注意事项 (Note)

For Research Use Only , Not for Diagnostic Use.