


概述 (Summary)

英文全称	Notch 1/Neurogenic locus notch homolog protein 1/Motch A/mT14
纯度 (Purity)	Greater than 95% as determined by reducing SDS-PAGE
内毒素 (Endotoxin level)	<1 EU/μg as determined by LAL test.
蛋白构建 (Construction)	Recombinant Mouse Notch 1 is produced by our Mammalian expression system and the target gene encoding Ala18-Gln526 is expressed with a 6His tag at the C-terminus.
Accession #	Q01705
蛋白标签 (Tag)	
表达宿主 (Host)	Human Cells
种属 (Species)	Mouse
预测分子量 (Predicted MW)	54.4 KDa
蛋白形态 (Form)	Lyophilized from a 0.2 μm filtered solution of PBS, pH 7.4.
储存缓冲液 (Buffer)	
运输方式 (Shipping)	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature listed below.
稳定性&储存 (Stability &Storage)	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
复溶 (Reconstitution)	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100μg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles. Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100μg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

电泳图 (SDS-PAGE image)

Product Name: Recombinant Mouse Notch 1 (C-6His)
Catalog #: PHM2018



背景 (Background)

分子别名 (Alternative Names)

Neurogenic locus notch homolog protein 1; Notch 1; Motch A; mT14

背景介绍 (References)

Mouse Notch1 is a 300 kDa type I transmembrane glycoprotein and it functions as a receptor for membrane-bound ligands Jagged1, Jagged2 and Delta1 to regulate cell-fate determination. Mouse Notch1 is synthesized as a 2531 amino acid (aa) precursor that contains an 18 aa signal sequence, a 1707 aa extracellular domain (ECD) with 36 EGFLike repeats and three Lin12/notch repeats, a 21 aa transmembrane segment and a 785 aa cytoplasmic domain that contains six ankyrin repeats, a glutamine-rich domain and a PEST sequence. Notch1 may play an essential role in postimplantation development, probably in some aspect of cell specification and/or differentiation and may be involved in mesoderm development, somite formation and neurogenesis.

注意事项 (Note)

For Research Use Only , Not for Diagnostic Use.