

产品名称: Recombinant Mouse SHH (C-6His)
产品货号: PHM2194




概述 (Summary)

英文全称	SHH/Sonic Hedgehog
纯度 (Purity)	Greater than 95% as determined by reducing SDS-PAGE
内毒素 (Endotoxin level)	<1 EU/μg as determined by LAL test.
蛋白构建 (Construction)	Recombinant Mouse Sonic Hedgehog is produced by our Mammalian expression system and the target gene encoding Cys25/xadGly198 is expressed with a 6His tag at the C-terminus.
Accession #	Q62226
表达宿主 (Host)	Human Cells
种属 (Species)	Mouse
预测分子量 (Predicted MW)	20.4 KDa
制剂 (Form)	Lyophilized from a 0.2 μm filtered solution of 20mM Citrate, 8% Trehalose, 4% Mannitol, 0.02% Tween80 (w/v), pH 5.5.
运输方式 (Shipping)	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature listed below.
稳定性&储存 (Stability &Storage)	Lyophilized protein should be stored at ≤ -20°C, stable for one year after receipt. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at ≤ -20°C for 3 months.
复溶 (Reconstitution)	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100μg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles. Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100μg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

电泳图 (SDS-PAGE image)

产品名称: Recombinant Mouse SHH (C-6His)
产品货号: PHM2194



背景 (Background)

分子别名 (Alternative Names)

Sonic Hedgehog Protein; SHH; HHG-1; SHH

背景介绍 (References)

Mouse Sonic Hedgehog Homolog (SHH) belongs to a three-protein family called Hedgehog. The other two family members are Indian Hedgehog (IHH) and Desert Hedgehog (DHH). Hedgehog proteins are key signaling molecules in embryonic development. SHH is expressed in various embryonic tissues and plays critical roles in regulating the patterning of many systems, such as limbs and brain. SHH also plays an important role in adult, including the division of adult stem cells and the development of certain cancers and other diseases. Mouse Shh is synthesized as a 437 aa precursor that contains a 24 aa signal sequence and a 413 aa mature region. The mature region is autocatalytically processed into a nonglycosylated, 20 kDa, 174 aa N/xadterminal fragment (Shh/xadN), and a catalytic/xadprocessing, glycosylated, 34 kDa, 239 aa C/xadterminal fragment. The 20 kDa Shh/xadN fragment is the core of the active hedgehog molecule. Mouse Shh/xadN is 99%, 98%, and 100% aa identical to human, rat and gerbil Shh/xadN, respectively.

注意事项 (Note)

For research use only .