


概述 (Summary)

英文全称	AXL/Tyrosine-protein kinase receptor UFO/AXL oncogene/UFO
纯度 (Purity)	Greater than 95% as determined by reducing SDS-PAGE
内毒素 (Endotoxin level)	<1 EU/μg as determined by LAL test.
蛋白构建 (Construction)	Recombinant Mouse Tyrosine-protein kinase receptor UFO is produced by our Mammalian expression system and the target gene encoding Ala19-Pro443 is expressed with a human IgG1 Fc tag at the C-terminus.
Accession #	Q00993
蛋白标签 (Tag)	
表达宿主 (Host)	Human Cells
种属 (Species)	Mouse
预测分子量 (Predicted MW)	73.3 KDa
蛋白形态 (Form)	Lyophilized from a 0.2 μm filtered solution of PBS, pH 7.4.
储存缓冲液 (Buffer)	
运输方式 (Shipping)	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature listed below.
稳定性&储存 (Stability &Storage)	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
复溶 (Reconstitution)	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100μg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles. Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100μg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

电泳图 (SDS-PAGE image)

Product Name: Recombinant Mouse AXL (C-Fc)
Catalog #: PHM2260



背景 (Background)

分子别名 (Alternative Names)

Tyrosine-protein kinase receptor UFO; AXL oncogene; UFO

背景介绍 (References)

Axl, also known as Ufo and Ark, is a widely expressed 140 kDa glycoprotein in the TAM receptor tyrosine kinase family. Axl binds the vitamin K-dependent protein Gas6 which triggers tyrosine autophosphorylation of the Axl cytoplasmic domain. Axl functions in dampening the immune response, regulating cytokine secretion, clearing apoptotic cells and debris, and maintaining cell survival. Axl is highly expressed in solid cancers and promotes in vivo tumorigenesis and tumor cell invasiveness. It also functions as a cellular entry receptor for Gas6-opsonized lentiviruses. Axl contributes to cell survival, migration, invasion, metastasis and chemosensitivity justify further investigation of Axl as novel therapeutic targets in cancer. The receptor tyrosine kinase AXL is thought to play a role in metastasis. The soluble AXL receptor as a therapeutic candidate agent for treatment of metastatic ovarian cancer. GAS6/AXL targeting as an effective strategy for inhibition of metastatic tumor progression in vivo.

注意事项 (Note)

For Research Use Only , Not for Diagnostic Use.