


Product Name: Recombinant Macaca nemestrina PCSK9 (C-6His)
Catalog #: PHV1855

概述 (Summary)

英文全称	PCSK9/Proprotein Convertase 9
纯度 (Purity)	Greater than 95% as determined by reducing SDS-PAGE
内毒素 (Endotoxin level)	<1 EU/μg as determined by LAL test.
蛋白构建 (Construction)	Recombinant Macaca Nemestrina Proprotein Convertase Subtilisin/Kexin Type 9/PCSK9 is produced by our Mammalian expression system and the target gene encoding Gln31-Gln152&Ser153-Gln692 is expressed with a 6His tag at the C-terminus.
Accession #	A8T662
蛋白标签 (Tag)	
表达宿主 (Host)	Human Cells
种属 (Species)	Macaca nemestrina
预测分子量 (Predicted MW)	1459 KDa
蛋白形态 (Form)	Supplied as a 0.2 μm filtered solution of 20mM PB, 150mM NaCl, 0.1M Arg, 0.1M Glu, 20% glycerol, 0.01% tween20, 5% Trehalose, pH 6.0.
储存缓冲液 (Buffer)	
运输方式 (Shipping)	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
稳定性&储存 (Stability &Storage)	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
复溶 (Reconstitution)	

电泳图 (SDS-PAGE image)



背景 (Background)

Product Name: Recombinant Macaca nemestrina PCSK9 (C-6His)
Catalog #: PHV1855

分子别名 (Alternative Names)

Proprotein Convertase Subtilisin/Kexin Type 9; Proprotein Convertase 9; PC9; Subtilisin/Kexin-Like Protease PC9; PCSK9

背景介绍 (References)

Proprotein Convertase Subtilisin/Kexin Type 9 (PCSK9) is a secretory subtilase belonging to the proteinase K subfamily. PCSK9 is synthesized as a soluble zymogen that undergoes autocatalytic intramolecular processing in the ER, the pro domain and mature chain secrete together through noncovalent interactions. PCSK9 binds with low-density lipoprotein receptor (LDLR) and plays a major regulatory role in cholesterol homeostasis. PCSK9 also plays a role in the neural development.

注意事项 (Note)

For Research Use Only, Not for Diagnostic Use.