


概述 (Summary)

英文全称	GITR Ligand/TNFSF18
纯度 (Purity)	Greater than 95% as determined by reducing SDS-PAGE
内毒素 (Endotoxin level)	<1 EU/μg as determined by LAL test.
蛋白构建 (Construction)	Recombinant Cynomolgus TNF Superfamily Member 18 is produced by our Mammalian expression system and the target gene encoding Glu74-Ser199 is expressed with a 6His tag at the C-terminus.
Accession #	A0A2K5UCD9
蛋白标签 (Tag)	
表达宿主 (Host)	Human Cells
种属 (Species)	Cynomolgus
预测分子量 (Predicted MW)	15.3 KDa
蛋白形态 (Form)	Lyophilized from a 0.2 μm filtered solution of PBS, pH 7.4.
储存缓冲液 (Buffer)	
运输方式 (Shipping)	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature listed below.
稳定性&储存 (Stability &Storage)	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
复溶 (Reconstitution)	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100μg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles. Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100μg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

电泳图 (SDS-PAGE image)

Product Name: Recombinant Cynomolgus GITR Ligand (C-6His)
Catalog #: PHV2078



背景 (Background)

分子别名 (Alternative Names)

Tumor necrosis factor ligand superfamily member 18; TNFSF18; Activation-inducible TNF-related ligand; Glucocorticoid-induced TNF-related ligand; AITRL; GITRL; TL6

背景介绍 (References)

TNFSF18 is a single-pass type II membrane protein. It is expressed at high levels in the small intestine, ovary, testis, kidney and endothelial cells. TNFSF18 cytokine binds to TNFRSF18/AITR/GITR. It regulates T-cell responses, and functions as costimulator and lower the threshold for T-cell activation and T-cell proliferation. It is Important for interactions between activated T-lymphocytes and endothelial cells and Promotes leukocyte adhesion to endothelial cells. TNFSF18 mediates activation of NF-kappa-B. As Triggers increased phosphorylation of STAT1 and up-regulates expression of VCAM1 and ICAM1. It also regulates migration of monocytes from the splenic reservoir to sites of inflammation.

注意事项 (Note)

For Research Use Only , Not for Diagnostic Use.